




What would you expect from a Pioneer in HbA1c?



A Complete Solution for HbA1c Testing Needs

As the number of diabetic patients increases, so does the requirement for HbA1c testing. You need a company that understands every aspect of HbA1c testing needs.



For HbA1c – Think Tosoh

TOSOH BIOSCIENCE



G8 HPLC ANALYZER
Gold Standard Method for HbA1c Measurement

G7 HPLC ANALYZER
HbA1c and β -Thalassemia in One Instrument



AIA-360
Compact and Affordable System for Flexible Testing

AIA-600 II
Full-Featured Benchtop Performance

HPLC HbA1c

Gold standard ion-exchange HPLC

1.6 minutes analysis time (G8)

Primary tube sampling for various tube sizes

Direct determination of stable HbA1c

Less than 2% CVs

Simple touch screen operation

Compact footprint (G8) – 21”h x 20”w x 19”d

Calibrators are directly traceable to the IFCC reference standard

G7 and G8 are NGSP certified


Low maintenance

No interference - up to 10% HbF

No interference from variants HbS, HbC and HbD

90 sample loader

290 sample loader and line automation models available



AIA HbA1c

Results are equivalent to Tosoh's gold standard HPLC

10 minutes incubation time

90 day calibration stability (helps minimize costs and reagent waste)

Direct determination of stable HbA1c (no separate measurement of Hb and HbA1c)

Less than 2% CVs (within assay)

HbA1c testing in 4 simple steps

Compact footprint (360) – 21”h x 16”w x 16”d

Calibrators are directly traceable to the IFCC reference standard

AIA-360 and AIA-600II are NGSP certified

Low maintenance

No interference - up to 10% HbF

No interference from variants HbC and HbD

Does not require a dedicated system for HbA1c

Uses Tosoh's well known unit dose test cup technology

A Complete Solution For HbA1c Testing Needs

Tosoh Provides Solutions




for your Growing HbA1c Testing Needs

Testing Volumes



Correlation



TOSOH BIOSCIENCE

Tosoh Bioscience, Inc.
 6000 Shoreline Court, Suite 101
 South San Francisco, CA 94080
 Tel: (800) 248-6764
 Fax: (650) 615-0415
www.tosohbioscience.us

Tosoh and AIA are registered trademarks of Tosoh Corporation. ©2011